

## Safety Data Sheet

### 1. PRODUCT AND COMPANY IDENTIFICATION

#### 1.1 Product identifiers

Product Name: **MicroSnap Enrichment Broth™**  
Product Number: MS1-CEC-BROTH-2ML  
Brand: Hygienea

#### 1.2 Recommended use and restrictions on use

Identified use: Enrichment Broth for Coliform and *E.coli* Test  
Restriction on use: Use per manufacturer's instructions only

#### 1.3 Details of the supplier of the safety data sheet

Company: Hygienea USA  
941 Avenida Acaso  
Camarillo, CA 93012  
USA  
Telephone: 805-388-8007  
Fax: 805-388-5531  
E-mail: info@hygienea.com

#### 1.4 Emergency telephone number

Emergency Phone: 1-888-494-4362

### 2. HAZARDS IDENTIFICATION

#### 2.1 Hazard Classification:

Not classified as hazardous according to OSHA Hazard Communication Standard 29 CFR 1910.1200, Regulation EC No. 1272:2008 [CLP/GHS] and WHIMS 2015 (HPR) Hazardous Products Regulation SOR 2015-17.

#### 2.2 Label elements: Not applicable

Pictograms: Not applicable

Hazard codes: Not applicable

#### 2.3 Hazards not otherwise clarified: None

### 3. COMPOSITION/INFORMATION ON INGREDIENTS

| Ingredient                                                                                                                                                                                | CAS #   | % by Wt |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------|---------|
| Non hazardous components in Broth                                                                                                                                                         | Mixture | 3%      |
| <b>FORM:</b> Aqueous nutrient media contained in plastic tube. When used according to manufacturer's instructions; the liquid stable reagents will not come in contact with the operator. |         |         |

### 4. FIRST AID MEASURES

#### 4.1 Description of first aid measures

#### 4.2 Inhalation:

If breathed in, move person into fresh air. If not breathing give artificial respiration and seek medical advice.

#### 4.3 Skin contact:

Wash off with soap and plenty of water. If you are concerned get medical advice.

#### 4.4 Eye contact:

Flush eyes with large amounts of water. Remove contact lenses if easy to do. Continue rinsing. If symptoms persist, get medical attention.

#### 4.5 If swallowed

Rinse mouth with water. If you feel unwell, get medical attention. Never give anything by mouth to an unconscious person

#### 4.6 Most important symptoms and effects, both acute and delayed

None anticipated.

---

**4. FIRST AID MEASURES (Continued)**

- 4.7 Indication of any immediate medical attention and special treatment needed**  
Not anticipated.
- 

**5. FIRE-FIGHTING MEASURES**

**5.1 Extinguishing media**

**Suitable extinguishing media**

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

**5.2 Special hazards arising from the substance or mixture**

None inherent in the products. Burning plastic will create carbon and nitrogen oxides.

**5.3 Advice for firefighters**

Wear self-contained breathing apparatus for firefighting if necessary.

**5.4 Further information**

No data available

---

**6. ACCIDENTAL RELEASE MEASURES**

**6.1 Personal precautions protective equipment and emergency procedure**

Ensure adequate ventilation. For personal protection see section 8.

**6.2 Environmental precautions**

Avoid release to the environment. Prevent further leakage or spillage if safe to do so. Do not let product enter drains.

**6.3 Methods and materials for containment and cleaning up**

Contain spill, collect using commercially available inorganic absorbent material and place in container for disposal according to local regulations (see section 13).

---

**7. HANDLING AND STORAGE**

**7.1 Precautions for safe handling**

Avoid contact with skin and eyes. Do not eat, drink or smoke when using this product. Wash hands after handling.

**7.2 Conditions for safe storage, including any incompatibilities**

Keep refrigerated between 2-8°C and keep in foil bag to protect from light.

**7.3 Specific end use(s)**

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

---

**8. EXPOSURE CONTROLS/PERSONAL PROTECTION**

**8.1 Control parameters**

None required.

**8.2 Exposure controls**

**Appropriate engineering controls**

Handle in accordance with good industrial hygiene and safety practice. Wash hands after use.

**Personal protective equipment Eye/face protection**

No eye protection required, safety goggles may be used for added protection.

**Skin/hand protection**

No protective gloves are required, but may be used for added protection. Wash hands after use.

**Respiratory protection**

Respiratory protection is not required for intended use.

**Control of environmental exposure**

Do not let product enter drains.

---

**9. PHYSICAL AND CHEMICAL PROPERTIES**

**9.1 Information on basic physical and chemical properties**

- |                   |                                                                                           |
|-------------------|-------------------------------------------------------------------------------------------|
| a) Appearance     | Form: Liquid in bulb of device, plastic tube with extrusion and solution.<br>Color: Clear |
| b) Odor           | No data available                                                                         |
| c) Odor Threshold | Not applicable                                                                            |

---

**9. PHYSICAL AND CHEMICAL PROPERTIES (Continued)**

|    |                                                    |                             |
|----|----------------------------------------------------|-----------------------------|
| d) | pH                                                 | Neutral                     |
| e) | Melting point/<br>freezing point                   | No data available           |
| f) | Initial boiling point<br>and boiling range         | No data available           |
| g) | Flash point                                        | No data available           |
| h) | Evaporation rate                                   | No data available           |
| i) | Flammability<br>(solid, gas)                       | Not applicable              |
| j) | Upper/lower<br>flammability or<br>explosive limits | No data available           |
| k) | Vapour pressure                                    | No data available           |
| l) | Vapour density                                     | No data available           |
| m) | Specific Gravity                                   | No data available           |
| n) | Water solubility                                   | Readily miscible with water |
| o) | Partition coefficient:<br>n-octanol/water          | No data available           |
| p) | Auto-ignition<br>temperature                       | Not applicable              |
| q) | Decomposition<br>temperature                       | No data available           |
| r) | Viscosity                                          | No data available           |
| s) | Explosive properties                               | No data available           |
| t) | Oxidizing properties                               | No data available           |

---

**10. STABILITY AND REACTIVITY**

**10.1 Reactivity**

Product is considered non reactive under normal intended use.

**10.2 Chemical stability**

Stable under recommended storage conditions.

**10.2 Possibility of hazardous reactions**

No data available

**10.4 Conditions to avoid**

Heat and direct sunlight.

**10.5 Incompatible materials**

No data available.

**10.6 Hazardous decomposition products**

Carbon and Nitrogen oxides from burning plastic.

---

**11. TOXICOLOGICAL INFORMATION**

**Skin corrosion/irritation**

No data available

**Serious eye damage/eye irritation**

No data available

**Serious eye damage/eye irritation**

No data available

**Respiratory or skin sensitisation**

No data available

**Germ cell mutagenicity**

No data available

---

---

**11. TOXICOLOGICAL INFORMATION (Continued)**

**Carcinogenicity**

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

**Reproductive toxicity**

No data available

**Specific target organ toxicity - single exposure**

No data available

**Specific target organ toxicity - repeated exposure**

No data available

**Aspiration hazard**

No data available

**Additional Information**

RTECS: Not available

---

**12. ECOLOGICAL INFORMATION**

**12.1 Toxicity**

No data available

**12.2 Persistence and degradability**

No data available

**12.3 Bioaccumulative potential**

No data available

**12.4 Mobility in soil**

No data available

**12.5 Results of PBT and vPvB assessment**

PBT/vPvB assessment not available as chemical safety assessment not required/not conducted

---

**13. DISPOSAL CONSIDERATIONS**

**Product**

Observe local, State and Federal environmental regulations.

**Contaminated packaging**

Dispose of as unused product.

---

**14. TRANSPORT INFORMATION**

**DOT (US)**

Not dangerous goods

**IMDG**

Not dangerous goods

**IATA**

Not dangerous goods

No special precautions needed for transport either within or outside premises.

---

**15. REGULATORY INFORMATION**

This safety datasheet complies with the requirements of OSHA Hazard Communication Standard 29 CFR 1910.1200, Regulation EC No. 1272:2008 [CLP/GHS] and WHIMIS 2015 (HPR) Hazardous Products Regulation SOR 2015-17.

**Safety, health and environmental regulations/legislation specific for the substance or mixture**

No data available

**SARA 302 Components**

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 302.

**SARA 313 Components**

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

**SARA 311/312 Components**

No SARA Hazard

**Massachusetts Right To Know Components**

No components are subject to the Massachusetts Right to Know Act.

---

---

**15. REGULATORY INFORMATION (Continued)**

**Pennsylvania Right To Know Components**

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by Pennsylvania Right To Know.

**New Jersey Right To Know Components**

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by New Jersey Right To Know.

**California Prop. 65 Components**

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

---

**16. OTHER INFORMATION**

The product is non hazardous according to OSHA Hazard Communication Standard 29 CFR 1910.1200, EC No. 1272:2008 [CLP/GSH] and WHIMS 2015 HPR Hazardous Products Regulation SOR 2015-17, and therefore the product is not considered hazardous.

The information on this SDS is based on information on SDS for ingredients provided by supplier.

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. Hygiena and its affiliates shall not be held liable for any damage resulting from handling or from contact with the above product

**Preparation Information**

Hygiena USA  
Product Safety  
1-888-494-4362

Issued Date: 8/9/16 Revision: C

---

**End of Safety Data Sheet**